P-096

The new method of the colony measurement in Ames test by using Digital Microscope Culturing System

Kentaro Tanaka¹, Mitsuyoshi Miyashita², Hiroyuki Ogawa², Hiroyuki Murakami¹, Sayuri Yamada¹, Rieko Kubo¹, Osamu Tajima¹

¹Kirin Holdings Company, Ltd. Central Laboratories for frontier Technology, Center for Food Safety Science; ²MicroBio Corporation

As a recent practice, the Digital Microscope Culturing System (DMCS), that automates the official analytical method, is widely used for quality control testing and food safety checks at many food-processing plants. Because the DMCS has the ability to measure the rate of growth of colonies on each dish and only detects colonies formed by live microorganisms, the system can measure colonies more precisely compared with conventional colony counters. Furthermore, because of this ability, DMCS automatically counts colonies even if there is foodstuff debris on the dish being tested. In this study, we investigated whether or not the DMCS is suitable to automate the enumeration in the Ames test by evaluating microorganisms TA100, TA98, TA1535, TA1537 and WP2uvrA. As a result, it was clearly observed that the DMCS detected and counted colonies even in the area adjacent to the edges of the dish. Furthermore, a good correlation was shown between the manual colony counting and automated colony counting processes in both the samples with foodstuff debris and those without. Therefore, we reached the conclusion that the DMCS was well suited to the Ames test and this particular method reduced the enumeration time and the labor-intensive work of manual colony counting at the same time.

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