How to use R2A Agar Media (for 1ml sample)

MicroBio Corporation

1. General Description


This prepared plated agar media is designed to culture microorganism that normally inhabits in water where nutrient for growth is scarce. It takes a sample of 1ml, to avoid the heat damage that might be caused if poured-plate method is used. This media also has a characteristic of keeping moisture for prolonged incubation.

2. Detection Procedure (spread-plate for 1ml sample)

2-1 Inoculate a media evenly with 1ml sample and leave it for 15 to 30 minutes.

2-2 Incubate the plates at 20 to 25 (or 30 to 35) degrees C temperature for 4 to 7 days.

When Biomatic DMCS is used, fully automated rapid detection and precise colony counts can be achieved.


3. Preservation


For preservation, keep and store the media in dark place at room temperature.

4. Detection Examples (Biomatic DMCS S-12 Data)

4-1 1ml Sample Spread-plate: *Pseudomonas fluorescens* (NBRC15842) on PPM024 media [DMCS S-12 Data (22 hour-detection) at 25 degrees C Incubation]

4-2 1ml Sample Spread-plate: *Methylobacterium extorquens* (NBRC15911) on PPM024 media [DMCS S-12 Data (80 hour-detection) at 25 degrees C Incubation]

4-3 Filter Method

Example) [DMCS S-12 Data at 25 degrees C Incubation]

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